ThermoFisher SCIENTIFIC

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The screenshot shows the Proteome Discoverer software interface. The 'Study Definition' window displays a list of samples (F1-F4) with IDs and names. The 'Analysis' window shows a workflow with 'Consensus Step' and 'Processing Step'.

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Mix 1:1

Trypsin digest

K containing peptide R containing peptide

SILAC
(Stable Isotope Labeling by Amino Acids in Cell Culture)

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


40

ThermoFisher
SCIENTIFIC

	SILAC	iTRAQ/TMT	LFQ	DIA
Quantify with	MS1 Precursor	MS2 Reporter	MS1 Precursor	MS2 Fragment
Coverage (1shot)	++	+++	++++	+++++
Coverage (Extend)	+++++ (Fraction)	+++++ (Fraction)		
Accuracy	+++++	++* (Ratio Compression)	+++	++++
Precision	++++	+++++	+++	++++
Throughput	+++ (Up to 3plex)	+++++ (Up to 16plex)	++ (1 by 1)	++ (1 by 1)
Easy to start	+++ (Cost)	+++ (Cost)	+++++	++++ (Processing tool)
Notes	Cell line only.	Compatible with any samples (tissue, fluid etc.) *with mobility or MS3 method, accuracy will be improved to +++ or ++++	Most simple method same as standard shotgun proteomics	PTM is not supported in library prediction

42

ThermoFisher
SCIENTIFIC

Proteome Discoverer software : Quan Proteomics

The screenshot displays the Proteome Discoverer interface with a table of protein identifications. Two red callout boxes highlight specific features:

- Ratios based on median values from replicates:** Points to the 'Median (Median)' column in the table.
- p-values can be calculated due to replicate measurements:** Points to the 'p-Value' column, which contains values like 1.76e-1.

Checked	Protein FDR	Master	Gene	Accession	Protein Name	Median (Median)	p-Value	Abundance	...
1	High	✓	ADP1	P27188	Argininosuccinate synthase [OS=Saccharomyces cerevisiae (strain ATCC 22016 / DSM 15204 / EC 4.1.3.15)]	3.223	1.76e-1	69.1	...
2	High	✓	TDH1	P03260	Glyoxaldehyde 3-phosphate dehydrogenase 1 [OS=Saccharo...	0.857		72.2	...
3	High	✓	GDV2	P49095	Glyoxal dehydrogenase (scapaboylating), mitochondrial [OS=...	1.454		88.6	...

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The collage illustrates various data visualization methods available in the software:

- PCA:** A scatter plot showing the separation of different experimental conditions.
- Volcano plots:** A plot of log2 fold change versus -log10 p-value, highlighting significant differentially expressed proteins.
- Pathway maps:** Network diagrams showing the biological context of identified proteins, such as mitochondrial membranes and organelle envelopes.
- Heat maps:** A heatmap with dendrogram showing relative protein abundance across multiple samples.
- Enrichment analysis:** A table showing the enrichment of specific gene ontology terms.

Accession	Description	Size	Overlap	Expected	Fold Enrichment	p-Value	-log2 (p-Value)
GO:0019888	organelle inner membrane	47	24	8.3	2.9	0.000105	23.26
GO:000743	mitochondrial inner membrane	47	24	8.1	2.9	0.000105	23.26
GO:0019866	mitochondrial membrane	68	29	12	2.4	0.000105	20.71

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Peptide ID	TMT SPSMS3	LFQ	Phosphorylation	Glycoproteomics	Crosslinking	TMT MS2
SILAC	Spectral Counting	Lipidomics	IBAQ	Intact Profiling	Middle-Down	DIA
Top Down	Denovo Sequencing	δ -carboxylation	Neu Code	Metabolic Labeling (SL)	Oxidation	AIM
PRM	Other PTMs	Combinatorial Isotopic Mass Tags	δ -Carboxylation	δ -glutamylation	δ -sulfonylation	Penylation
Ubiquitylation	ADP Ribosylation	SUMOylation	SUMOylation	SUMOylation	SUMOylation	SLIM
Metabolomics	Glycation	Glycation	Glycation	Glycation	Glycation	SLIM
ITRAQ MS2	WiSIM-DIA	Lasso peptides	SIM	Acylation	Acylation	DipyrO
Dimethyl Labeling	S-nitrosylation	Disulfide Structures	Limited Proteolysis	Pegylation	Disulfide Mapping	Dileu Labeling
Glycomics	iTRAQ SPSMS3	Peptide Epimer Analysis	iodoTMT	Cysteine Residue	S-Sulfination	Thiol Modifications
Acetylation	Peptide Mapping	S-Palmitoylation	Hydroxylation	Thiol Modifications	Thiol Modifications	Hydroxylation

- Peptide ID
- Quan Proteomics
- Phosphorylation

45

Proteomics : Deep : Phosphopeptide Enrichment

DEEP

Fe-NTA

3,378 Proteins

12,465 Phosphopeptides

10,436 Quantified

46